




Mayer's Hematoxylin Solution

REF CS-0001-20

 40 (20 ml)

For use in chromogenic *in situ* hybridization procedures



In vitro diagnostic medical device
according to EU directive 98/79/EC

1. Intended use

The Mayer's Hematoxylin Solution (CS1) is intended to be used for counterstaining steps in chromogenic *in situ* hybridization (CISH) procedures. The Mayer's Hematoxylin Solution is intended to be used in combination with *ZytoDot* or *ZytoFast* kits (T-1151, C-3003, C-3005 or C-3018).

Interpretation of the results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

2. Clinical relevance

Refer to the instruction for use of the respective *ZytoDot* or *ZytoFast* kit.

3. Test principal

The chromogenic *in situ* hybridization (CISH) technique allows for the detection and visualization of specific nucleic acid sequences in cell preparations. Digoxigenin-labeled DNA or RNA fragments, so called CISH probes, and their complementary target DNA or RNA strands in the preparations are co-denatured and subsequently allowed to anneal during hybridization. Duplex formation of the digoxigenin-labeled probe and the target DNA strand is indirectly detected by using a digoxigenin-binding antibody which in turn is detected by a secondary polymerized enzyme-conjugated antibody. The enzymatic reaction of a chromogenic substrate leads to the formation of a color precipitate. After counterstaining the DNA with a nuclear dye, color precipitates are visualized using a light microscope.

4. Reagents provided

Mayer's Hematoxylin Solution (CS-0001-20):

- Quantity: 20 ml
- Amount of tests: sufficient for 40 tests

5. Materials required but not provided

- *ZytoDot* or *ZytoFast* kit (T-1151, C-3003, C-3005 or C-3018)
- Microscope slides, positively charged
- Water bath (55°C, 80°C, boiling)
- Hot plate or hybridizer
- Humidity chamber + hybridization oven or hybridizer
- Adjustable pipettes (10 μ l, 1000 μ l)
- Staining jars or baths
- Timer
- Calibrated thermometer
- Ethanol or reagent alcohol
- Xylene
- Deionized or distilled water
- Coverslips (22 mm x 22 mm, 24 mm x 32 mm)
- Rubber cement, e.g., Fixogum Rubber Cement (Prod. No. E-4005-50/-125) or similar
- Adequately maintained light microscope

6. Storage and handling

The component must be stored at 2-8°C. If these storage conditions are followed, the component will function, without loss of performance, at least until the expiry date printed on the label.

7. Warnings and precautions

- Read the instruction for use prior to use!
- Do not use the reagents after the expiry date has been reached!
- Do not reuse reagents!
- This product contains substances (in low concentrations and volumes) that are harmful to health and potentially infectious. Avoid any direct contact with the reagents. Take appropriate protective measures (use disposable gloves, protective glasses, and lab garments)!
- If reagents come into contact with skin, rinse skin immediately with copious quantities of water!
- A material safety data sheet is available on request for the professional user.
- The specimens must not be allowed to dry during the hybridization and washing steps!
- Avoid any cross-contamination and micro-bacterial contamination of the reagents!

For further information concerning this point please refer to the instruction for use of the respective *ZytoDot* or *ZytoFast* kit.

8. Limitations

- For *in vitro* diagnostic use.
- For professional use only
- The clinical interpretation of any positive staining, or its absence, must be done within the context of clinical history, morphology, other histopathological criteria as well as other diagnostic tests. It is the responsibility of a qualified pathologist to be familiar with the CISH probes, reagents, diagnostic panels, and methods used to produce the stained preparation. Staining must be performed in a certified, licensed laboratory under the supervision of a pathologist who is responsible for reviewing the stained slides and assuring the adequacy of positive and negative controls.
- Specimen staining, especially signal intensity and background staining, is dependent on the handling and processing of the specimen prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning, or contamination with other specimens or fluids may produce artefacts or false results. Inconsistent results may result from variations in fixation and embedding methods, as well as from inherent irregularities within the specimen.
- The performance was validated using the procedures described in the instruction for use of the respective *ZytoDot* or *ZytoFast* kit. Modifications to these procedures might alter the performance and have to be validated by the user.

9. Interfering substances

Refer to the instruction for use of the respective *ZytoDot* or *ZytoFast* kit.

10. Preparation of specimens

Refer to the instruction for use of the respective *ZytoDot* or *ZytoFast* kit.

11. Preparatory treatment of the device

The device is ready-to-use. No reconstitution, mixing, or dilution is required.

12. Assay procedure

For detailed information on how to perform chromogenic *in situ* hybridizations with *ZytoDot* and *ZytoFast* products, including counterstaining steps with the Mayer's Hematoxylin Solution (CS1), please refer to the instruction for use of the respective *ZytoDot* or *ZytoFast* kit.

13. Interpretation of results

Refer to the instructions for use of the respective *ZytoDot* or *ZytoFast* kit.

14. Recommended quality control procedures

Refer to the instructions for use of the respective *ZytoDot* or *ZytoFast* kit.

15. Disposal

The disposal of reagents must be carried out in accordance with local regulations.

16. Troubleshooting

Any deviation from the operating instructions can lead to inferior staining results or to no staining at all. Please refer to the instructions for use of the respective *ZytoDot* or *ZytoFast* kit for further information.

17. Literature

- Hopman AHN, et al. (1997) *Histochem Cell Biol* **108**: 291-8.
- Isola J, Tanner M (2004) *Methods Mol Med* **97**: 133-44.
- Shipley J (2006) *J Pathol* **210**: 1-2.
- Speel EJ, et al. (1994) *J Histochem Cytochem* **42**: 1299-307.
- Tsukamoto T, et al. (1991) *Int J Dev Biol* **35**: 25-32.
- Wilkinson DG: *In Situ Hybridization, A Practical Approach*, Oxford University Press (1992) ISBN 0 19 963327 4.

Trademarks:

ZytoVision[®], *ZytoDo*[®] and *ZytoFast*[®] are trademarks of *ZytoVision GmbH*.



ZytoVision GmbH
Fischkai 1
27572 Bremerhaven/ Germany
Phone: +49 471 4832-300
Fax: +49 471 4832-509
www.zytovision.com
Email: info@zytovision.com

Our experts are available to answer your questions.
Please contact helptech@zytovision.com