

Wash Buffer SSC

REF WB-0001-560 ∇ 8 (560 ml)

For use in *in situ* hybridization procedures




IVD

In vitro diagnostic medical device

according to EU directive 98/79/EC



Wash Buffer SSC

REF WB-0001-560  8 (560 ml)

For use in *in situ* hybridization procedures

CE

IVD

In vitro diagnostic medical device

according to EU directive 98/79/EC

1. Scope of Application

This product is designed for *in vitro* diagnostic use (according to EU directive 98/79/EC). Interpretation of results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of patient by a qualified pathologist!

The Wash Buffer SSC (WB1) is designed to be used for washing steps in *in situ* hybridization (ISH) protocols of ZytoVision.

2. Safety Precautions and Disposal

- ✓ Read the operating instructions prior to use!
- ✓ Do not use the reagents after the expiry date has been reached!
- ✓ Avoid any cross-contamination and micro-bacterial contamination of the reagents!
- ✓ Avoid any direct contact with the reagents. Take appropriate protective measures (use disposable gloves, protective glasses, and lab garments)!
- ✓ If reagents come into contact with skin, rinse skin immediately with copious quantities of water!
- ✓ Never pipet solutions with your mouth!
- ✓ The disposal of reagents must be carried out in accordance with local regulations!
- ✓ A material safety data sheet is available on request for the professional user!

3. The Wash Buffer SSC

The following components are included:

Code	Component	Quantity	Container
WB1	<u>Wash Buffer SSC</u>	560 ml	Screw-cap bottle (large)
	Instruction manual	1	

Component (**WB1**) is sufficient for 8 staining jars of 70 ml each.

4. Storage and Shelf Life

The component must be stored at 2...8°C. If these storage conditions are followed, the component will function, without loss of performance, at least until the expiry date printed on the label

5. Instructions

Wash Buffer SSC (WB1) is a ready-to-use solution for washing steps in ISH applications, following the ISH protocols of ZytoVision.

Our experts are available to answer your questions.


As of: June 15, 2018 (4.6)

Trademarks:

ZytoVision® is a trademark of ZytoVision GmbH.



ZytoVision GmbH · Fischkai 1
D - 27572 Bremerhaven · Germany
Phone: +49 (0)471/4832 - 300
Fax: +49 (0)471/4832 - 509
www.zytovision.com
info@zytovision.com



Your local distributor